



## O2k-Info

Mitochondrial Physiology Network 19.17: 1-2 (2014)  
Updates: [http://wiki.orooboros.at/index.php/MiPNet19.17\\_Series\\_G](http://wiki.orooboros.at/index.php/MiPNet19.17_Series_G)

©2014 OROBOROS®  
Version 01: 2015-04-09

# Oxygraph-2k Series G: innovations

Fleischmann S, Fasching M, Gradl P, Schwaninger H, Gnaiger E

**OROBOROS INSTRUMENTS Corp.**  
high-resolution respirometry  
Schöpfstr. 18, A-6020 Innsbruck, Austria  
[instruments@orooboros.at](mailto:instruments@orooboros.at); [www.orooboros.at](http://www.orooboros.at)



The current version of the O2k-Core (OROBOROS Oxygraph-2k Core) - the experimental system for basic HRR - is the Oxygraph-2k Series G. The O2k-Series G sets a new standard in merging high-end technology and Corporate Social Responsibility.

### The added features include:

- Increased energy efficiency thanks to the new O2k-electronics design: lower running costs (24 W compared to 41 W in Series E-F and 136 W in Series D at 37 °C) following our general concept for sustainable technology. The max. power input is 120 W.
- O2k net weight is reduced to 13.45 kg in comparison to 14.2 kg in Series E-F and 21 kg in Series D.
- [DatLab 6](#) is the new software version for the Oxygraph-2k Series G.
- Ambient temperature: An internal thermistor is implemented for continuous monitoring of ambient room temperature, as an additional element of quality control. This sets a new standard in high-resolution respirometry, (i) at the level of automatic feedback-control elements in the firmware and (ii) by allowing the experimenter to evaluate any possible influence of external temperature fluctuations on signal stability.
- External temperature port (TEMP ext; Pt1000): An external thermistor can be connected, e.g. for monitoring the temperature in the aqueous medium of the O2k-Chamber. This is a relevant objective when designing respiratory measurements in response to dynamic temperature changes. The additional temperature channels can be continuously recorded by DatLab 6.
- Faster connection time to DatLab.

