


http://wiki.oroboros.at/index.php/O2k-Publications: Lymphocyte
High-resolution respirometry and cryopreservation of lymphocytes

## Mitochondrial respiration in peripheral blood mononuclear cells correlates with depressive subsymptoms and severity of major depression



Karabatsiakis A, Böck C, Salinas-Manrique J, Kolassa S, Calzia E, Dietrich DE, Kolassa IT

## High-resolution respirometry of cryopreserved lymphocytes as a diagnostic tool for major depression (MD)



**Figure 1.** Boxplots of group-wise comparison of respiration in PBMCs from MD patients (N = 19) and control subjects (N = 19) characterized by mitochondrial oxygen flux ( $I_{02}$ ) per cell. Circles and triangles represent statistical outliers of the respective groups. Asterisk (\*) indicates significant group differences on an alpha level of 0.05 (one-tailed Mann–Whitney U test or Student's independent t-test). MD, major depression; PBMCs, peripheral blood mononuclear cells; SRC, spare respiratory capacity

The significant decrease of ET-capacity and excess E-R capacity factor indicates that cells of acutely depressed patients have a reduced capacity of mitochondrial respiration that contribute to the pathophysiology and/or aetiology of depressive symptoms

Reference: Karabatsiakis A, Böck C, Salinas-Manrique J, Kolassa S, Calzia E, Dietrich DE, Kolassa IT (2014) Mitochondrial respiration in peripheral blood mononuclear cells correlates with depressive subsymptoms and severity of major depression. Transl Psychiatry 4:e397.

Figures and texts slightly modified based on the recommendations of the COST Action MitoEAGLE CA15203. Doi:10.26124/mitofit:190001.v4

## O2k-brief communicated by LF Garcia-Souza Oroboros Instruments

