

OROBOROS O2k-Info



Mitochondrial Physiology Network 07.07(04):1 (2017)

Version 04: 2017-06-29 ©2017 OROBOROS

Updates: http://wiki.oroboros.at/index.php/MiPNet07.07_O2k-NetworkInvitation

Dear colleague,

It is our honour to invite you to the

O2k-Network

as a Reference Laboratory for O2k high-resolution respirometry.

International collaboration, support and development are the keystones of this progressive international network that comprises more than 350 laboratories. Through joining this network, customer support by the OROBOROS-Team and our research laboratory in Innsbruck can be complemented by direct user-to-user contacts.

The **O2k-Network** offers:

- Direct information on the methodological details of Oxygraph-2k applications used by other O2k-Network labs.
- The prospect of combining the expertise of leading research laboratories in the field of O2k high-resolution respirometry.

As an **O2k-Network Reference Laboratory** for O2k high-resolution respirometry:

- Your lab is listed on the OROBOROS website¹ with information on your team and areas of research.
- Your O2k-Publications regarding the OROBOROS O2k are listed on your O2k-Network page.
- You receive advance information on O2k-Workshops and specific OROBOROS Events.
- We keep you alert on international projects within and beyond the O2k-Network: [MitoEAGLE](#), [K-Regio MitoFit](#)

We are looking forward to your positive response. With your cooperation, the O2k-Network will become an even bigger success!

Erich and the OROBOROS-Team, instruments@oroboros.at

¹<http://wiki.oroboros.at/index.php/O2k-Network>